

# Filtering variants in a clinical setting

MSc in Genomic Medicine

Lucy Crooks

26/1/2016

## SDGS Gene panels

- 4 Panels of 300 kb (CTD, Hereditary cancers, NGD, IEM)
- 70-100 genes
- 5 million paired-end reads up to 1000X coverage
- 96 samples per run
- Confidence that no variants are missed

# SDGS variant analysis

Alignment

BWA    Picard    GATK

Variant calling

GATK HaplotypeCaller

Base  
quality $\geq$ 10  
Mapping  
quality $\geq$ 20

Quality filtering

Remove variants  
with QUAL $<$ 30

Identify key variant

VCFtools (filter to target)

Remove polymorphisms (in-house list)

SNPEff

ANNOVAR

Sanger sequence gaps with  
coverage  $<$  30X

# Pipeline Validation

- 100% detection rate for 300 variants found by Sanger sequencing
  - 31 patients
  - 28 genes
  - 3 panels