





Sequencing & Analysing Cancer Genomes


Challenges of Variant Calling in Cancer

Matthew Parker, Ph. D
Lead Bioinformatician
Sheffield Diagnostic Genetics Service

- Clinical
 - Turn around times
 - Appropriate samples
 - Variant prioritisation
 - "Action-ability"
- Analysis
 - Quality of DNA
 - Tumor Content
 - Clonality
 - BMT




How do you think variants discovered in a cancer genome could be used?




Mike Stratton

Genomics England Cancer Programme



Vectors by Vecteezy.com

- Sample quality affects downstream sequence quality
- Problems:
 - Difficult to get enough DNA sometimes
 - PCR May be required → Not ideal for WGS
 - Necrotic tissue – poor quality DNA
 - Storage method affects quality



NATURE COMMUNICATIONS | ARTICLE [OPEN](#)




Systematic pan-cancer analysis of tumour purity

Dvir Aran, Marina Sirota & Atul J. Butte

[Affiliations](#) | [Contributions](#) | [Corresponding author](#)

Nature Communications 6, Article number: 8971 | doi:10.1038/ncomms9971



Received 02 September 2015 | Accepted 19 October 2015 | Published 04 December 2015



Adenoid cystic carcinoma

Lung Adenocarcinoma

Liver Hepatocellular Carcinoma





REVIEW

doi:10.1038/nature10762

Clonal evolution in cancer

Mel Greaves¹ & Carlo C. Maley²




REVIEW

doi:10.1038/nature10762

Clonal evolution in cancer

Mel Greaves¹ & Carlo C. Maley²

- Imagine 80% tumor purity.
- Coverage of 10x
- 8/10 reads on average are from the tumor
- Imagine 3 clones, one at 50%, two at 25%
- Het Mutations – only one read in the two smaller clones
- Only 2 reads supporting variant in largest clone




Coverage is King!

- Evidence of tumor material in the normal sample is sometimes evident
- Tumor cell infiltration of surrounding tissues/blood/saliva
- Some callers can't handle this and the SNV is excluded (for example Illumina!!)

| GeneName | SJQuality | Chr | Pos | Tumor MAF | Normal MAF | Total Normal Coverage |
|----------|-----------|-------|-----------|-----------|------------|-----------------------|
| HSPA12A | SJHQ | chr10 | 118464762 | 0.29 | 0.05 | 21 |
| ZNF365 | SJHQ | chr10 | 64415123 | 0.31 | 0.00 | 16 |
| AKR1C1 | SJHLQ | chr10 | 5005711 | 0.32 | 0.04 | 25 |
| GRAMD1B | SJLQ | chr11 | 123448167 | 0.13 | 0.00 | 14 |
| GIT2 | SJHQ | chr12 | 110389058 | 0.29 | 0.06 | 35 |
| KCNH5 | SJHQ | chr14 | 63174997 | 0.34 | 0.07 | 29 |
| KLC1 | SJLQ | chr14 | 104139491 | 0.09 | 0.00 | 21 |
| HMG2P46 | SJHQ | chr15 | 45808097 | 0.33 | 0.06 | 18 |
| NR2E3 | SJHQ | chr15 | 72109938 | 0.28 | 0.05 | 21 |
| DSG1 | SJHQ | chr18 | 28936281 | 0.30 | 0.04 | 25 |
| NUDC | SJHQ | chr1 | 27250594 | 0.34 | 0.03 | 31 |
| ISM1 | SJHLQ | chr20 | 13251341 | 0.30 | 0.09 | 35 |
| TRMU | SJHQ | chr22 | 46733794 | 0.29 | 0.04 | 24 |
| CADPS | SJHQ | chr3 | 62459955 | 0.42 | 0.06 | 18 |
| MTTP | SJHQ | chr4 | 100485294 | 0.32 | 0.09 | 33 |
| ZNF718 | SJHQ | chr4 | 154730 | 0.31 | 0.06 | 31 |
| SETP8 | SJHQ | chr5 | 132096615 | 0.33 | 0.05 | 20 |
| LARP1 | SJLQ | chr5 | 154195506 | 0.10 | 0.00 | 23 |
| PCDHA8 | SJLQ | chr5 | 140223016 | 0.20 | 0.00 | 8 |
| PTPRK | SJHQ | chr6 | 128290476 | 0.31 | 0.09 | 34 |
| ARHGEF10 | SJHQ | chr8 | 1772199 | 0.40 | 0.00 | 18 |
| OR5C1 | SJHQ | chr9 | 125551552 | 0.23 | 0.04 | 28 |
| GPR143 | SJHQ | chrX | 9733666 | 0.34 | 0.00 | 15 |
| NHS | SJHQ | chrX | 17743982 | 0.24 | 0.04 | 26 |

Evidence of Contamination?




The patients blood is now the blood of the donor – using as a normal sample for somatic calls – not possible

- Because of Purity Heterogeneity – need a relatively high coverage Tumor sample > 30x
- Normal sample can be 20-30x
- Example Projects:
 - 1000G
 - PCGP
 - TCGA
- GEL: Tumor >75x, Normal 35x


- Mapping is the same – usually BWA
- Variant calling is done by algorithms specifically designed for cancer
 - Not A diploid genomes
 - Tumor admixture in the normal
 - Subclonal variants
 - Copy number variants, copy-neutral LOH or major ploidy changes
- Trade-off between sensitivity (detecting variants) and specificity (rate of false positives)
- Mutect is a good example of a cancer variant caller

- Variants for the cancer (somatic) are called with the germline sample in mind – at the same time.
- You cannot just do a simple subtraction
- Generally algorithms do some simple stats:
 - Probability of the variant being real given the reads from cancer sample and the reads from the germline sample
 - Usually involves post processing of data – filtering etc
 - Mutations found at high % in the population for example are usually filtered out

- SNVs in cancer allow assessment of:
 - Sample-sample contamination
 - Tumour & Normal provenance – same patient?
 - Numbers/types of variants




Contamination Example



- Aligners usually have a difficult time with indels
- Usually happen in more repetitive regions of DNA
- Left alignment


```
CCTCGGCCCTCCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTg
CCTCGGCCCTCCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTC
CCTCGGCCCTCCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTC
CCTCGGCCCTCCTCCTCCTCCTCCTCGGC---CCTCCTCCTC
CCTCGGC---CCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTC
CCTCGGCCCTCCTCCTCCTCCTCCTCGGC---CCTCCTCCTC
CCTCGGCCCTCCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTC
CCTCGGCCCTCCTCCTCCTCCTCCTCGGCCCTCCTCCTCCTC
```


- 3 sources of information can be used to detect SVs:
 - Soft-clipped reads
 - Discordant read pairs
 - Unmapped Reads




- Example: CREST:
<http://www.stjuderresearch.org/site/lab/zhang>

- Read depth information is leveraged



CONSERTING: integrating copy-number analysis with structural-variation detection

Xiang Chen^{1,2}, Pankaj Gupta^{1,2}, Jianmin Wang²⁻⁴, Joy Nakitandwe^{2,5}, Kathryn Roberts⁵, James D Dalton⁵, Matthew Parker^{1,2}, Samir Patel⁵, Linda Holmfeldt⁵, Debbie Payne⁵, John Easton^{2,6}, Jing Ma^{2,5}, Michael Rusch^{1,2}, Gang Wu^{1,2}, Aman Patel^{1,2}, Suzanne J Baker^{2,7}, Michael A Dyer^{2,7}, Sheila Shurtleff^{2,5}, Stephen Espy³, Stanley Pounds⁸, James R Downing^{2,5}, David W Ellison^{2,5}, Charles G Mullighan^{2,5} & Jinghui Zhang^{1,2}



Parker et al. *Genome Biology* 2012, **13**:R113
<http://genomebiology.com/2012/13/12/R113>




RESEARCH


Open Access

Assessing telomeric DNA content in pediatric cancers using whole-genome sequencing data

Matthew Parker¹, Xiang Chen¹, Armita Bahrami², James Dalton², Michael Rusch¹, Gang Wu¹, John Easton³, Nai-Kong Cheung⁴, Michael Dyer⁵, Elaine R Mardis^{6,7}, Richard K Wilson^{6,7}, Charles Mullighan², Richard Gilbertson⁵, Suzanne J Baker⁵, Gerard Zambetti⁸, David W Ellison², James R Downing² and Jinghui Zhang^{1*}, for the Pediatric Cancer Genome Project

- Deregulation of telomere maintenance is a hallmark of cancer.
- Most tumors find a way to maintain telomere length above a critical minimum
- Can be measured in WGS data





*silent mutations removed



Whole Genome Analysis from 100,000 Genomes Project Cancer Programme

Release notes for version 1.1 (composite WGA):


- WGAs are currently being returned from the Initiation Implementation Phase of the Cancer Main Programme on participants with solid tumours for whom paired fresh frozen and germline samples were submitted and passed quality assurance checks.
- It is envisaged that the initial review of these WGAs will be carried out within the GMC accredited Molecular Pathology/Genetics Laboratory therefore these analyses will be returned to GMC nominated individuals in these laboratories.
- WGAs are currently identified via two IDs; the Genomics England Participant ID and the Lab Sample ID for the fresh frozen tumour sample in line with NHS working practice. Provision of additional identifiers is under evaluation.
- The WGA is released as a paired set of formats:
 - a **preliminary analysis** containing small variants in Domain 1 (77 genes containing variants annotated as potentially actionable by Genome Oncology) and Domain 2 (590 genes listed in the Cancer Gene Census).
 - a **supplementary analysis** containing the content of the preliminary analysis and additional small variants (occurring outside Domains 1 & 2), copy number variants

The screenshot shows the My Cancer Genome website in a web browser. The browser's address bar displays the URL <https://www.mycancergenome.org>. The website's header includes the My Cancer Genome logo and a search bar. Below the header, there are navigation tabs for Home, DIRECT, and About Us. The main content area is divided into two columns. The left column contains two sections: 'Find a Cancer Mutation' and 'Find Clinical Trials'. The 'Find a Cancer Mutation' section has three dropdown menus for Disease (required), Gene (optional), and Variant (optional), all with 'Select Disease First' as the default selection, and a 'GO' button. The 'Find Clinical Trials' section has two text input fields for Disease (optional) and Gene (optional), and a 'GO' button. The right column contains three sections: 'Learn About My Cancer Genome' with a right-pointing arrow, 'Support My Cancer Genome' with a 'GIVE' button and a 'More...' link, and 'Molecular Medicine' with a list of links including Pathways, Immunotherapy in Cancer, Overview of Targeted Therapies for Cancer, List of Anticancer Agents, Types of Molecular Tumor Testing, and Circulating Tumor DNA. Below these is a 'Connect With Us' section with a list of links including Newsletter, Newsletter Sign-Up, MCG Mobile App, and Take Our Survey. The footer of the website contains a navigation bar with links to MCG Home, About Us, Acknowledgments, Give, Site Map, Legal, and APIs and Licensing. At the bottom, it states 'My Cancer Genome is managed by the [Vanderbilt-Ingram Cancer Center](#) Copyright © 2010 - 2017 MY CANCER GENOME' and 'Powered by VANDERBILT'.

Genetically Informed Can x

Secure | <https://www.mycancergenome.org>

Apps SCH Inbox (2,791) MED676 sch-sdgs BioInf HUB NGS-Server Tmux Projects Explain SAM Flags Swagger UI bcbio GATK | GATK | Tool D Catalog OLAT Other bookmarks

 **MY CANCER GENOME**
GENETICALLY INFORMED CANCER MEDICINE

Search My Cancer Genome

Home DIRECT About Us

Find a Cancer Mutation

Disease (required):

Gene (optional):

Variant (optional):

Find Clinical Trials


Lists trials by Disease or Gene for all trials registered within the U.S. National Cancer Institute's [list of NCI-supported clinical trials](#).

Disease (optional):

Gene (optional):

Learn About My Cancer Genome ▶

Support My Cancer Genome


 Help create new tools and resources

[More...](#)

Molecular Medicine

- ▶ Pathways
- ▶ Immunotherapy in Cancer
- ▶ Overview of Targeted Therapies for Cancer
- ▶ List of Anticancer Agents
- ▶ Types of Molecular Tumor Testing
- ▶ Circulating Tumor DNA


Connect With Us



- ▶ Newsletter
- ▶ Newsletter Sign-Up
- ▶ MCG Mobile App
- ▶ Take Our Survey

MCG Home | About Us | Acknowledgments | Give | Site Map | Legal | APIs and Licensing

My Cancer Genome is managed by the [Vanderbilt-Ingram Cancer Center](#) Copyright © 2010 - 2017 MY CANCER GENOME

Powered by  VANDERBILT

Whole Genome Analysis

100,000 Genomes Project Cancer Programme

Preliminary analysis of somatic small non-synonymous variants v1.1



Participant information

| Participant name | D.O.B | Gender | NHS number | Laboratory sample ID | GeL participant ID | GMC | Sample date | Date analysis issued |
|------------------|-------|--------|------------|----------------------|--------------------|-----|-------------|----------------------|
| XX | | | | | | | | |

Tumour information

| Tumour type | Tumour subtype | ICD10 code | Sample type | Reported tumour content | Tumour sample cross-contamination |
|-------------|----------------|------------|-------------|-------------------------|-----------------------------------|
| Colorectal | adenocarcinoma | N/A | FF | Medium 40-60% | PASS |

Domain 1 variants

Variants in a virtual panel of potentially actionable genes*. Actionable genes are defined as genes in which small variants (SNVs and indels <50bp) have reported therapeutic, prognostic or clinical trial associations**, as defined by the GenomOncology Knowledge Management System. Where known, the "variant-level actionability" category and applicable tumour type are indicated. For other variants in these genes, their impact on gene function has not yet been characterised and therefore their actionability status is unclear. This means:

- (i) local evaluation will be required for listed variants which are not yet characterised (i.e. "variant-level actionability" is denoted N/A)
- (ii) even if well characterised as actionable for some tumour types, the listed variants may not be actionable in the participant's specific tumour type

*Current potentially actionable genes for solid tumours: 77 genes, listed at [Actionable genes in solid tumour v1.1](#) document

**Links are provided to clinical trials within the United Kingdom which are both actively recruiting participants or closed to recruitment.

Right Tumour, Reported Actionability.

| Gene | Gene-level actionability | GRCh38 coordinates ref/alt allele | Transcript | cDNA and protein change | Predicted consequences | Population germline allele frequency (1KG) | VAF | Alt allele/ total read depth | COSMIC ID | Variant-level actionability | Gene mode of action |
|------|--|-----------------------------------|-----------------|---------------------------|------------------------|--|------|------------------------------|--|---|-----------------------------|
| ALK | Therapeutic (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (solid neoplasm); Trial (solid neoplasm); Trial (solid neoplasm) | 2:29220810 G>A | ENST00000389048 | c.3541C>T p. (Arg1181Cys) | missense_variant | N/A | 0.15 | 20/130 | N/A | Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (solid neoplasm); Trial (solid neoplasm); Trial (solid neoplasm) | oncogene |
| KRAS | Therapeutic (colorectal ca); Therapeutic (NSC lung ca); Trial (colorectal ca); Trial (colorectal ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (NSC lung ca); Trial (solid neoplasm); Trial (solid neoplasm); Trial (solid neoplasm); Trial (solid neoplasm); Trial (glioma); Trial (MPNST); Trial (melanoma); Trial (neuroblastoma); Trial (rhabdoid tu); Trial (rhabdomyosarcoma); Trial (schwannoma); Trial (sarcoma-ST) | 12:25225628 C>T | ENST00000256078 | c.436G>A p. (Ala146Thr) | missense_variant | N/A | 0.3 | 35/118 | COSM19404 COSM1165198 | Therapeutic (colorectal ca); Therapeutic (NSC lung ca); Trial (colorectal ca); Trial (NSC lung ca); Trial (colorectal ca); Trial (glioma); Trial (MPNST); Trial (melanoma); Trial (neuroblastoma); Trial (NSC lung ca); Trial (rhabdoid tu); Trial (rhabdomyosarcoma); Trial (schwannoma); Trial (sarcoma-ST); Trial (solid neoplasm); Trial (solid neoplasm); Trial (solid neoplasm) | oncogene |
| TP53 | Trial (ovarian ca) | 17:7673796 C>A | ENST00000269305 | c.824G>T p. (Cys275Phe) | missense_variant | N/A | 0.49 | 40/82 | COSM10701 COSM99932 COSM3723938 COSM1637959 | Trial (ovarian ca) | oncogene, tumour suppressor |

Domain 2 variants

Variants in a virtual panel of cancer-related genes***. Cancer-related genes are defined as genes in which any variants have been causally implicated in cancer, as defined by the Cancer Gene Census (Wellcome Trust Sanger Institute)

***Current cancer-related genes: 590 genes, listed at [Cancer census genes v1.1](http://cancer.sanger.ac.uk/census/) document

Genes causally implicated in cancer: Cancer Gene Census

<http://cancer.sanger.ac.uk/census/>

| Gene | GRCh38 coordinates ref/alt allele | Transcript | cDNA and protein change | Predicted consequences | Population germline allele frequency (1KG) | VAF | Alt allele/total read depth | COSMIC ID | Gene mode of action |
|--------|-----------------------------------|-----------------|-----------------------------------|------------------------|--|------|-----------------------------|--|-----------------------------|
| APC | 5:112839499 TGCAA>T | ENST00000508376 | c.3906_3909delGCAA p.(Leu1302>fs) | frameshift_variant | N/A | 0.13 | 14/110 | N/A | tumour suppressor |
| CREBBP | 16:3758051 T>TA | ENST00000262367 | c.3370-4dupT | splice_region_variant | N/A | 0.2 | 12/61 | N/A | oncogene, tumour suppressor |
| FAT1 | 4:186620732 C>T | ENST00000441802 | c.5854G>A p.(Val1952Ile) | missense_variant | N/A | 0.28 | 28/99 | COSM1054196 COSM1054194 | tumour suppressor |
| FIP1L1 | 4:53453080 CAG>C | ENST00000337488 | c.1459_1460delAG p.(Arg483>fs) | frameshift_variant | N/A | 0.11 | 9/83 | COSM249696 COSM4435275 | N/A |
| IDH2 | 15:90088607 T>C | ENST00000330062 | c.514A>G p.(Arg172Gly) | missense_variant | N/A | 0.24 | 27/112 | COSM33731 | oncogene |

Domain 3 variants

Small variants in genes not included in domains 1 & 2. These are not included in this document but are accessible via the Supplementary Analysis.

Sequencing quality information

See online [Technical Information v1.1.main](#) document and/or LabKey QC portal for details and expected ranges of QC metrics

| Sample type | Mapped reads, % | Chimeric DNA fragments, % | Insert size median, bp | Genome-wide coverage mean, x | Unevenness of local genome coverage, x | COSMIC content with low coverage (<30x), % | Total somatic SNVs | Total somatic indels | Total somatic SVs |
|-------------|-----------------|---------------------------|------------------------|------------------------------|--|--|--------------------|----------------------|-------------------|
| Germline | 95.38 | 0.40 | 482.8 | 29.3 | 6.65 | N/A | N/A | N/A | N/A |
| Tumour | 95.68 | 0.37 | 447.4 | 84.6 | 13.45 | 1.15 | 31152 | 21184 | 262 |

Additional information

- The pathways for sample processing and data analysis are not yet accredited end-to-end for diagnostic use. Accordingly, any result intended for use in informing clinical management should be confirmed using a test accredited for clinical use.
- Sensitivity: the depth of WGS used in this analysis will detect 99% of SNVs with an allele frequency of ≥ 0.3 , 95% of SNVs with an allele frequency of ≥ 0.1 and 60% of indels with an allele frequency of ≥ 0.4 (estimate is based upon admixtures analysis of a highly accurate catalog of variants produced in the "platinum genomes" project). Consequently, variants with allelic frequencies below this level, or in areas of low coverage may not be detected. False negative results cannot be excluded.
- Somatic calls are filtered according to the quality and quantity of reads. Full details of the filters used in this analysis can be found in the [Technical Information v1.1.main](#).
- Variants present in the germline are subtracted to produce a list of somatic variants. Accordingly, variants detected in both the germline and the tumour will not be listed in this analysis.

- In this analysis MNVs (multiple nucleotide variants) are reported as multiple consecutive SNVs and therefore the protein change may require correction.
- Only variants with specific consequences (transcript ablation, splice acceptor variant, splice donor variant, stop gained, frameshift variant, stop lost, initiator codon variant, transcript amplification, inframe insertion, inframe deletion, missense variant, splice region variant, incomplete terminal codon variant) in canonical transcripts are reported. The complete list of canonical transcripts can be accessed at [List of canonical transcripts v1.1](#).
- A variant may have multiple entries in COSMIC database due to the use of different reference sequences. In these cases links to all COSMIC entries are provided.
- Structural variants (SVs) and copy number variants (CNVs) are not included in this analysis. These variants types are included in the Supplementary Analysis.
- For a full description of the methods used to produce these results and for further information regarding QC metrics please refer to the [Technical Information v1.1.main](#). All related documentation is available at [Genomics England Website](#).
- 'N/A' indicates that information is not available or not applicable.

Genomics England

Queen Mary University of London
Dawson Hall
Charterhouse Square
London
EC1M 6BQ

Sequencing Laboratory

Illumina Laboratory Services United Kingdom - Hinxton
The Ogilvie Building, Wellcome Trust Genome Campus
Hinxton Nr Saffron Walden
Essex
CB10 1DR

GEL Supplementary Cancer Report

- Utility of cancer genomes in the clinic not yet realised
- Need good samples with high tumour content
- Coverage is King!
- Variant calling is challenging
 - Heterogeneity
 - Tumour purity
 - Contamination
- Knowledgebases for prioritising variants increasingly important – “Actionability”